

PRODUCT INFORMATION

Clinical Center Standardization Committee



GLOVE UPDATE

July 2004 (2)

REGENT MEDICAL BIOGEL® REPORTS A SHORT-LIVED MANUFACTURING CHANGE TO THE STERILE SKINSENSE PI GLOVES.

THE 'GRIP RATE' AT THE FINGERTIPS HAS BEEN INCREASED FOR A LIMITED TIME.

CUSTOMERS MAY NOTICE A SUBTLE TACTILE SENSATION CHANGE.

THE GLOVE THICKNESS IS THE SAME (10 ML).

CUSTOMERS WILL BE NOTIFIED IN A FUTURE 'PRODUCT INFORMATION' WHEN THE 'GRIP RATE' IS BACK TO THE ORIGINAL SPECIFICATION.



MANAGEMENT TIPS FOR THE TEMPORARY INCREASED 'GRIP RATE':



DO NOT DOUBLE GLOVE WITH TWO SKINSENSE PI.



DOUBLE GLOVE WITH EITHER SKINSENSE UNIVERSAL AS FIRST GLOVE AND SKINSENSE PI AS SECOND GLOVE OR TWO SKINSENSE UNIVERSAL GLOVES.



USE THE SKINSENSE UNIVERSAL (ONLY 8 ML THICK) IF THE PI INCREASED GRIP RATE INTERFERES A LOT WITH TACTILE SENSATION.



UNDERSTAND THAT THIS INCONVENIENCE IS UNAVOIDABLE AND WILL BE SHORT-LIVED.



Skinsense™ PI (light blue, sterile, surgery, latex-free, smooth finish, powder-free, polyisoprene, 10 ml thick)

Item ID#: 02810 [Size 5.5]; 02811 [6.0]; 02812 [6.5]; 02813 [7.0]; 02814 [7.5]; 02815 [8.0]; 02816 [8.5]; 02817 [9.0]

Find in VSC: Glove, Surg, Latex Free, Biogel, Skinsense, [glove size], [item#]



Skinsense™ N Universal (dark blue, sterile, surgery, latex-free, smooth finish, powder-free, neoprene, 8 ml thick)

Item ID#: 02818 [Size 5.5]; 02819 [6.0]; 02820 [6.5]; 02821 [7.0]; 02822 [7.5]; 02823 [8.0]; 02824 [8.5]; 02825 [9.0]

Find in VSC: Glove, Surg, Latex Free, Biogel, Skinsense N Univ, [glove size], [item#]