

PRODUCT INFORMATION

Clinical Center Standardization Committee



UPDATE

BAC T/ ALERT BLOOD CULTURE BOTTLE & BLOOD CULTURE BOTTLE ADAPTER DECEMBER 2003 (01)

✧ If you have any of the BacT/ Alert bottles with the wide top (about one inch diameter) left in your inventory please use these bottles.

- Use with clear blood culture adapter ONLY.
- If you have no clear adapters then call micro 301-496-4433 to switch the bottles to the new slender style.
- The clear adapter has been discontinued and is no longer available from CHS or the Visual Supply Catalog.



✧ If you have the new BacT/ Alert bottles (slim top, about one half inch) you can use either blue adapter or the clear adapter.

- CHS provides the blue adapter on Visual Supply Catalog (adapter blood culture, Bac-T Alert, 00453)



either



or



✧ Micro supplies BacT/ Alert blood culture bottles. Call 301-496-4433 if you need blood culture bottles.

✧ Remember to send blood cultures in paired (aerobic and anaerobic) sets.