

**CLIA RESOURCE CENTER  
CLINICAL PATHOLOGY DEPARTMENT  
NATIONAL INSTITUTES OF HEALTH  
BLDG 10, RM 2C306  
BETHESDA, MD 20892**

**CLIA CHECKLIST**

**Quality Assurance**

- Lab director or designee reviews all results
- Inconsistent results are investigated and resolved with documentation in writing
- Testing personnel are documented on an annual basis to be competent in all tasks performed
- Review QC data record lot #'s expiration dates, open dates
- Review PT data and document
- Keep discontinued procedure for 2 yrs.
- Communicate changes to staff (staff meetings)
- Validation of in house tests-precision/accuracy: sensitivity/specificity
- Corrected reports-contact requester with corrected report alert limits report
- Review results for relevance to patient data (i.e. age, sex, and other pertinent clinical data)
- Action taken if discrepancy is discovered, monitor complaints received from users
- QA records are maintained for at least two years
- Keep discontinued procedures for two years

**Quality Control**

- Laboratory had appropriate environmental conditions (i.e. lab has defined condition in which testing may be performed, temperature, humidity, etc.)
- Procedure manuals

All tests procedures are written and complete including

- Test name
  - Principal of the test
  - Necessary equipment
  - Directions for specimen collection and handling
  - Directions for calibration or standardization
  - Directions for preparation of reagents, standards and controls
  - Step-by-step directions for performing the test
  - Quality control procedures and criteria defining unacceptable control results
  - Corrective actions when control criteria are exceeded
  - Reference ranges
  - Directions for calculation of results where appropriate
  - Notes, safety procedures
  - References
- 
- New procedures are reviewed and signed by director
  - Procedures are reviewed by the lab director or designee annually
  - All tests must be validated prior to reporting patient results. Precision, accuracy,

- sensitivity, and specificity of the test must be established.
- Reference ranges (if pertinent) must be established and reported with the results
  - A written QC policy includes:
    - Frequency of performing controls for each test/record lot #'s and date
    - At least two levels of control samples must be tested once a day of operation with each run (i.e. high /low or positive/negative)
    - The type of control
    - Acceptable limits for the control results
    - Results not reported if out of control
  - Corrective actions taken when results have exceeded acceptable limits have been reviewed
  - A policy to assure abiding by QC criteria is established
  - Instrumentation
    - Equipment is calibrated every 6 months and documentation is available for review
    - Equipment is maintained according to manufacturer's recommendations
    - Maintenance documents are reviewed
    - Temperature dependent instruments are checked and recorded each day of use
    - Acceptable ranges are defined for all temperature dependent equipment
    - Minimum /mid point/ maximum values required to verify range
    - Comparison of results of same specimen ran on two instruments of the same model, make, etc.
  - Reagents, dated, expiration date, etc. when open
  - Keep records two yrs.

### **Patient Test Management**

- Written policies and procedures are available for specimen collection, labeling, and preservation and handling including collection container, vol., transport criteria, labeling (safeguards for same first, last name)
- Written criteria are available for rejection of unacceptable specimens
- Test requisition includes:
  - Name of patient with identifier
  - Name and address of individual ordering test
  - Test to be performed
  - Date (time if critical) of specimen collection
  - Any additional information relevant and necessary to assure accurate and timely testing
  - Oral requests must have written follow up within 30 days
- Experimental data
  - Test name
  - Patient name or unique ID
  - Date received and/or tested
  - Test results
  - Indicate if unacceptable specimen
  - QC results and acceptable ranges
  - Indicate that QC results are acceptable

- Indicate who performed test
- Instrument printout of patient test results
- Test reports to physicians and or patients
  - Lab name and address of laboratory performing test
  - Patient name or unique ID
  - Test name
  - Test results
  - Normal ranges
  - Indicate if unacceptable specimen (i.e. condition and disposition of specimen)
  - Date of report
  - Special interpretation notes explanation
  - Disclaimer
    - "This test was developed and its performance characteristics determined by (insert name of lab). It has not been cleared or approved by the U.S. Food and Drug Administration.
- Test results secured confidential
- Easily retrieved
- If specimen is not tested--report includes disposition of specimen (repeat specimen requested)
- If incomplete information is received steps to be taken -documents attempts to receive needed information, and attempts to prevent reoccurrence
- Lab maintains record of reports for 2 years on a manner that can be readily retrieved
- Lab maintains record of information contributing to test result for 2 years on a manner that can be readily retrieved
- Referral of specimens for testing to a laboratory possessing a valid CLIA certificate
- Expected turn around time

### **Proficiency Testing**

- Testing is performed every six months using:
  - Commercially available analyte
  - Or-
  - Positive and negative samples from previous runs are performed every 6 months or high and low level sample for quantitative test
- Samples are blinded to the tester
- Results are held maintained and available for review for a minimum of 2 years
- Established target values for PT material
- Acceptable levels of variation are established and results meet acceptable or unacceptable criteria
- Document corrective action- notify testing personnel to prevent reoccurrence
- Review for effective corrective action- no repeat occurrence
- Review of proficiency testing: initial of Director/date

### **Personnel Qualifications**

- Education and experience of all personnel are on file in the human resources department or a copy of highest degree may be found in the laboratory standard operating procedure manual

- Personnel policies are available for review
- Testing personnel have had competency documented once a year; new employees at six months and annually thereafter
- Training of new personnel is documented

**Safety**

- Universal precautions are followed
- Training is documented annually
- Radiation Safety (RSB)
- Environmental monitors- room temp 65-75
- If particular requirements -state acceptable limits and how they are monitored