

# PRODUCT ALERT

Clinical Center Standardization Committee



MAY 2003 (3)

## GEMSTAR® EPIDURAL LINE PUMP SET DISCONTINUED

### SITUATION

Effective immediately the GemStar Epidural Line Pump Set—Item ID #02766—is discontinued.

The design configuration of the Pharmacy “PRN” IV bags is not compatible with Item ID #02766.

*“The discontinued tubing has an IV spike tip that cannot be used with the “PRN” IV bags.”*



### RECOMMENDED ACTION FOR OUR CUSTOMERS

➤ Do not use Item ID #02766.

➤ Use Item ID #02763  
GemStar Pump Set SL

UNIT OF ISSUE:

EACH

HOW TO FIND:

SET GEMSTAR ADMIN W/

INJECTOR



➤ Remember to identify tubing with ‘EPIDURAL LINE ONLY’ label.

➤ Place all unopened epidural line pump set—Item ID #02766—into soiled utility room. CHS will collect product during routine rounds;

➤ A summary of the Gemstar system infusion pump tubing (Table One) is attached.

➤ A new epidural line pump set design is under development. An announcement will be made when this is available.

# PRODUCT ALERT

Clinical Center Standardization Committee



MAY 2003 (3)  
TABLE ONE



## ABBOTT GEMSTAR® SYSTEM INFUSION PUMP TUBING SUMMARY

### OLD TUBING

TUBING NOW USED. USE WILL NOT CONTINUE.

#### ITEM ID# 00627

IV EXT ABBOTT PUMP

#### ITEM ID# 01452

ANESTHESIA PUMP SET

#### ITEM ID# 00416

PUMP SET W/ QUICK  
LOAD CARTRIDGE

#### ITEM ID # 02766

SET GEMSTAR LFSHLD  
PUMP 110"

### NEW TUBING

#### ITEM ID#

#### 02765

SET GEMSTAR PUMP

#### ITEM ID #

#### 02764

SET GEMSTAR LFSHLD  
PUMP 96"

#### ITEM ID #

#### 02763

SET GEMSTAR ADMIN W/  
INJECTOR

### TUBING NOW USED.

USE WILL CONTINUE.

#### ITEM ID # 00626

SET EXT ABBOTT PUMP  
IV EXT SL PROVIDER

#### ITEM ID # 00057

IV EXT FILTER 15"  
HIGH PRESSURE FILTER IV

#### ITEM ID # 00533

SET PROVIDER PUMP, 0.2 FILTER