

PRODUCT INFORMATION


Clinical Center Standardization Committee



Hemoccult Sensa Product Clarification August 2003 (1)

Three Hemoccult Sensa Product Are Available:

- Hemoccult Sensa Slide for Screening (Item ID 00122)
- Hemoccult Sensa Developer—Individual Size (Item ID 02598)
- Hemoccult II Sensa Patient Screening Kit—Dispak Plus (Item ID 02778)



Two of the three products are restricted to Hemoccult Point of Care Testing approved units:

Hemoccult Sensa Slide for Screening (Item ID 00122). This is ordered per box. Each box has 100 slides, two bottles of developer and 40 specimen applicators. This is the kit for Point of Care fecal occult bleeding testing.

Hemoccult Sensa Developer—Individual Size (Item ID 02598). This is ordered by individual vial. The kit for Point of Care fecal occult bleeding testing contains two bottles of developer, so extra developer may not be needed. Approved Point of Care units should order additional developer only as needed. It is important to maintain a limited supply on the unit for the purpose of temperature and to ensure that expired vials are not used.

Units approved for Hemoccult Point of Care Testing are 2 East, OP1, OP9, OP11 and OP12.

The third product may be requested by any PCU, including units approved for Hemoccult Point of Care Testing.

This is a patient take home kit/ envelope/ mailer containing instructions, three slides for fecal occult blood screening, and return envelope. Be sure to complete the return address to the Clinical Center Health Care Provider. Upon receipt of the mailer, the Clinical Center Health Care Provider enters medical order for test into computer and sends the mailer to DLM-Microbiology.

This is the Hemoccult II Sensa Patient Screening Kit—Dispak Plus (Item ID 02778).

Summary

		
Hemoccult Sensa Slide for Screening (Item ID 00122)	Hemoccult Sensa Developer—Individual Size (Item ID 02598)	Hemoccult II Sensa Patient Screening Kit—Dispak Plus (Item ID 02778)
Only for Hemoccult Point of Care Testing approved units.	Only for Hemoccult Point of Care Testing approved units.	May be requested by any PCU

NOTE: DLM-MICROBIOLOGY ACCEPTS ONLY THESE PRODUCTS FOR FECAL HEMOCCULT PROCESSING.