




# PRODUCT INFORMATION

## Clinical Center Standardization Committee



### PRODUCT REVIEW—STOPCOCKS

MAY 2004 (1)

-  An ORS stated that IV fluid (lipid infusion) was leaking via a stopcock.
-  Some stopcocks are Lipid Resistant.
-  Some stopcocks are Not Lipid Resistant.

---

#### LIPID RESISTANT STOPCOCKS




Item # 00807. Catalog # 42383-01. A 3-way Sorenson Intralock stopcock.

How to find in VSC: stopcock, 3-way, intralock, 00807.



Item # 00537. Catalog # 473545. A 4-way stopcock. Also nitroglycerin compatible. How to find in VSC: stopcock, 4-way, 4way, accel, 00537.



Item # 01043. Catalog # 456020. A 4-way stopcock. How to find in VSC: stopcock, 4-way, 4 way, 01043.

---

#### NOT LIPID RESISTANT STOPCOCKS



Item # 00049 Catalog # 041220001A. A double 4-way stopcock. How to find in VSC: stopcock, 4-way, 4 wy, double, 00049.



Item # 00074. Catalog # 70015004. A double Morse 3-way stopcock. How to find: stopcock, 3-way, double, 00074.




Item # 00232. Catalog # 70015003. A double Morse 3-way stopcock. How to find: stopcock, 3-way, 3 way w port, 00232.



Item # 02059. Catalog # 70015000. A 3-way stopcock. How to find in VSC: stopcock, 3 way, 3 way w/o port, 02059.



Item # 01120. Catalog # 42385-01. A 4-way Sorenson Intralock stopcock. How to find in VSC: stopcock, 4-way, 4 way, soren, 01120.



Item # 00324. Catalog # 2C6233. A 3-way Novex stopcock. How to find in VSC: stopcock, 3-way, novex, 3 way, 00324.

---