

PRODUCT UPDATE

Clinical Center Standardization Committee



CHANGED


BacT/ ALERT Bottle Shape Blood Culture Adapter Shape & Color December 2003 (1)

BacT/ ALERT Bottle Neck Shape Change:

From about one inch diameter to a half inch diameter.

CALL MICRO 301-496-4433 TO SWITCH OUT YOUR BOTTLES IF YOU HAVE BOTTLES WITH THE WIDER NECK.

Micro supplies BacT/ Alert bottles. BacT/ Alert bottles are not on Visual Supply Catalog or in CHS.




New bottle neck is slender.

Blood Culture Adapter Shape & Color Change Discontinued Item:

➔ About 2 inches tall. Clear color. **Only use with the one inch diameter bottles.**

➔ Do not use clear adapter with new half inch diameter bottles.



Replacement Item:

➔ About 2 ¼ inches tall. Blue color. **Only use with the new slender necked bottles.**

➔ Do not use blue adapter with wider neck bottles.

CHS Provides Adapter

Adapter Blood Culture
Item ID#: 00453
Unit of Issue: Each

How to find in VSC:
Adapter Blood Culture,
Bac-T Alert, 00453