

# PRODUCT UPDATE

Clinical Center Standardization Committee



April 2003 (1)

## AVAILABLE FROM CHS

### Chemotherapy Dispensing Pin

The BD Safety-Med chemotherapy dispensing pin has a 0.22 micron filter and an extra-sharp point for easy penetration.

**Unit of Issue:** Each

**Item ID#:** 02737

**How to find in MIS:**

P—Pin—Chemo

**How to find in VSC:**

Dispensing Pin, Chemosafety  
Dispensing Pin



**IMPORTANT Usage Tips:**

- Insert the pin **straight down without twisting** (this is due to the oval shape of the dispensing pin). DO NOT TWIST THE PIN.
- Failure to insert pin properly may cause leaking.

## ALERT UPDATE ON STOOL CONTAINER

RECENTLY WE REPORTED A DEFECT WITH OUR STOOL SPECIMEN CONTAINER. THE DEFECT WAS RELATED TO A SINGLE MANUFACTURED LOT.

WE HAVE REMOVED ALL DEFECTIVE PRODUCT. WE HAVE NEW PRODUCT IN STOCK THAT HAS BEEN TESTED TO ENSURE PRODUCT ADEQUACY.



CUSTOMERS NOW MAY RESUME USE OF THE STOOL SPECIMEN CONTAINER.

**Item ID#:** 00106

**How to find in MIS:**

C—Container—Stool

**How to find in VSC:**

Container Stool, Specimen Cup w/ Screw Lid