

PRODUCT INFORMATION

Clinical Center Standardization Committee




SAFETY ENHANCEMENT

NOVEMBER 2004 (1)

THE GLASS GRAY TOP 5 ML URINE TUBE IS AVAILABLE IN A PLASTIC RED/GRAY SPECKLED TOP 4 ML URINE TUBE.

CONTINUE TO USE THE GLASS GRAY TOP UNTIL SUPPLY IS GONE.

TECHNIQUE AND USE OF THE NEW PLASTIC TUBE IS THE SAME.



ITEM ID: 02434

UNIT OF ISSUE: EACH

HOW TO FIND IN VSC:

TUBE, URINE, PLASTIC, 5 ML, DRAW, GRAY, 02434